

ProSource® TF20

ENFit compatible, designed for tube feeding use



- Contains no sugars or artificial sweeteners, flavorings or colors.

FEATURES AND SPECIFICATIONS

Kcal/mL..... 0.89

Caloric Distribution (% of kcal)

Protein 100%
Carbohydrate 0%
Fat..... 0%

Protein Source:

Bovine Collagen

PDCAAS 100: Complete protein with 100% Dispensable and Indispensable amino acids.

HCPCS Code: B4155

Appropriate for these diets:

Lactose Free, Gluten Free, Diabetic,
Low Residue, Fluid Restricted

Shelf Stable

Do not refrigerate

Use immediately after opening.

Shelf Life: 18 months

ProSource TF20 Liquid Protein is thin as water and can be given directly down a feeding tube without mixing for easy administration and less risk of contamination.

No mixing or dilution required. Easy to dispense. ProSource TF20 Liquid Protein will not clog a clean well maintained feeding tube.

Not recommended to be added directly to enteral formula.

Not a sole source of nutrition.

A Medical Food

Use under medical supervision.

NOT FOR PARENTERAL USE

Ordering information:

Item Number	Case UPC	Billing Code	Packaging
11446 Unflavored	5-07-94688-11446-7	94688-0114-46	60/60 mL ENFit pouches

FOR THE DIETARY MANAGEMENT OF PROTEIN DEFICIENCY

ProSource® TF20 Liquid Protein is designed for tube feeding administration to meet the protein needs for critical illness, protein energy malnutrition and protein deficiency with good GI tolerance and greater patient application.

Recent studies indicate that critically ill patients require a minimum of 1.2-1.5 g/kg/day of protein and as much as 2+ g/kg/day of protein.^{41,42,43} Most enteral formulas available are able to supply adequate calories, but they are unable to meet the protein target without overfeeding. Subsequently, tube fed patients, often require a protein supplement to meet their individual protein needs, but not all protein supplements are appropriate for all patients.

Administration: Flush tube with 30 mL of water before and after administration.

For ENFit® Feeding Tubes: While holding pouch in an upright position, unscrew cap and insert into the ENFit connector. Invert pouch and gently squeeze to drain. When empty, disconnect and discard the pouch.

For Non-ENFit Feeding Tubes: Dispense into cup, draw up into syringe and administer. Appropriate for all enteral tubes. **Do not mix with tube feeding formula. Single use only. For enteral feeding only. Not a sole source of nutrition. Use under medical supervision.**

Ingredients: Water, Beef Collagen Hydrolyzate, Amino Acid Blend (L-Leucine, L-Phenylalanine, L-Lysine, L-Methionine, L-Isoleucine, L-Valine, L-Histidine, L-Threonine, L-Tryptophan), Natural Flavor.

Nutrition Facts

1 serving per container
Serving Size: 1 Pouch (60 mL)

Amount Per Serving

Calories 80

% Daily Values*

Total Fat 0g **0%**

Saturated Fat 0g **0%**

Trans Fat 0g

Cholesterol 0mg **0%**

Sodium 30mg **1%**

Potassium 4mg **0%**

Phosphorus 3mg **0%**

Total Carbohydrate 0g **0%**

Dietary Fiber 0g **0%**

Total Sugars 0g

Includes 0g Added Sugars **0%**

Protein 20g

Vit. D 0mcg 0% • Calcium 0mg 0%

Iron 0mg 0%

* The % Daily Value (DV) tells you how much a nutrient in a serving of food contributes to a daily diet. 2,000 calories a day is used for general nutrition advice.